

Development of a Multi-Wavelength Spatial Domain Cerebral Oximeter for Pediatrics.

L.A. Nelson MD, J.C. McCann, A.W. Loepke MD, C.D. Kurth MD

Institute of Pediatric Anesthesia, Cincinnati Children's Hospital Medical Center and the University of Cincinnati College of Medicine, Cincinnati, OH USA

Introduction: Detection of cerebral hypoxia-ischemia in infants remains problematic, as current monitors are impractical, insensitive, or nonspecific. Near infrared spectroscopy (NIRS) is capable of detecting cerebral hypoxia-ischemia at the bedside, although uncertainties persist about the accuracy of commercially available devices. Our study was designed to develop and test a new multi-wavelength spatial domain NIRS device in a piglet cerebral hypoxia-ischemia model.

Methods: The Institutional Review Committee approved the study. The prototype NIRS device consists of main unit connected to an optical probe housing photodiode detectors 2-4 cm from a 3 wavelength light emitting diode. An algorithm was derived from a Beer-Lambert spatial domain construct. Piglets were anesthetized with fentanyl (10µg/kg/hr infusion) and midazolam (0.1mg/kg/hr). Arterial and sagittal sinus saturation (SaO₂, SsO₂) was determined by co-oximetry. The algorithm was tested in 8 piglets subjected randomly to varying conditions of hypoventilation, hypoxia (low inspired oxygen concentration), and ischemia (carotid occlusion). Arterial blood gases, SaO₂ and SsO₂ were recorded at each condition. Mean cerebral saturation (SmO₂) was $SmO_2 = 0.85SsO_2 + 0.15SaO_2$ based on a 15:85 arterial:venous brain blood ratio. NIRS cerebral saturation (ScO₂) was evaluated relative to SmO₂ by bias, precision, slope, intercept, and correlation coefficient.

Results: We found a linear relationship ($ScO_2 = 0.72 SmO_2 + 14$) with good correlation ($r^2 = 0.84$). The bias and precision were 0% and 10%, respectively (Figure).

Discussion: This prototype NIRS device uses inexpensive, robust, clinically friendly technology similar to pulse-oximetry. The bias and precision is similar to other clinical devices such as pulse-oximetry, although accuracy remains intact over a wide saturation range. The multiwavelength spatial domain approach appears to improve accuracy. Clinical studies are underway.




Figure: Comparison of NIRS and arterial/venous oxygen saturation ratio

ScO₂ refers to cerebral oxygen saturation as measured with NIRS, SmO₂ refers to mean cerebral oxygen saturation, a 15:85 ratio of arterial and sagittal sinus venous oxygen saturation